

---

# **Epics Applications using PyEpics**

***Release 1.0***

**Matthew Newville**

March 23, 2012



# CONTENTS

<b>1</b>	<b>AreaDetector Display</b>	<b>3</b>
<b>2</b>	<b>Strip Chart</b>	<b>5</b>
<b>3</b>	<b>Epics Instruments</b>	<b>7</b>
<b>4</b>	<b>Sample Stage</b>	<b>9</b>
<b>5</b>	<b>Motor Setup</b>	<b>11</b>
<b>6</b>	<b>Ion Chamber</b>	<b>13</b>
<b>7</b>	<b>XRF Collector</b>	<b>15</b>
7.1	Overview . . . . .	15
7.2	Downloading and Installation . . . . .	15
7.3	Area Detector Display . . . . .	16
7.4	Strip Chart Display . . . . .	17
7.5	Using Epics Instruments . . . . .	19
7.6	Using Epics Motor Setup . . . . .	22
7.7	XRF Collector . . . . .	24
7.8	Ion Chamber . . . . .	24



PyEpics Apps contains several Epics Applications written in python, using the pyepics module (see <http://pyepics.github.com/pyepics/>). Many of these are GUI Application for interacting with Epics devices through Channel Access. The programs here are meant to be useful as end-user applications, or at least as examples showing how one can build complex applications with PyEpics. Many of the applications here rely on wxPython, and some also rely on other 3rd party modules (such as Image and SQLAlchemy).

The list of applications should be expanding, but currently include:



# **AREADETECTOR DISPLAY**

An application to control and view images from an Epics Area Detector. The controls available in this viewer are minimal, but you can change mode, exposure time and frame rate, and start and stop the Area Detector. The image can be manipulated by zooming, rotating, and so on.





# STRIP CHART

A simple “live plot” of the recent history of a set of PV values, similar to the common Epics StripTool. Time ranges and ranges for Y values can be changed, and data can be saved to plain text files.



# EPICS INSTRUMENTS

This application helps you organize PVs, by grouping them into instruments. Each instrument is a loose collection of PVs, but can have a set of named positions that you can tell it. That is, you can save the current set of PV values by giving it a simple name. After that, you can recall the saved values of each instrument, putting all or some of the PVs back to the saved values. The Epics Instrument application organizes instruments with tabbed windows, so that you can have a compact view of many instruments, saving and restoring positions as you wish.



# SAMPLE STAGE

A GSECARS-specific GUI for moving a set of motors for a sample stage, grabbing microscope images from a webcam, and saving named positions.



# MOTOR SETUP

An application for setting up and saving the configuration of Epics Motors. For each opened motor, a full setup screen is shown in a tabbed notebook display. A paragraph for a Motors.template file can be saved for each motor, or copied to the system clipboard. In addition, you can save and read “known motor types” to a database, which holds most of the motor parameters.

The Local version uses a local sqlite database to store and read the known motor types, while the GSE version (which works only at GSECARS, of course) connects to a mysql server on the GSECARS network to store and read the known motor types.





# ION CHAMBER

This non-GUI application is synchrotron-beamline specific. It reads several settings for the photo-current of an ion chamber and calculates absorbed and transmitted flux in photons/sec, and writes these back to PVs (an associated .db file and medm .adl file are provided). The script runs in a loop, updating the flux values continuously.



# XRF COLLECTOR

This non-GUI application interacts with a small epics database to save data from a multi-element fluorescence detector. The script runs as a separate process, watching PVs and saving data from the detector on demand. Associated .db file and medm .adl file are provided.

## 7.1 Overview

An *Epics Instrument* is simply a grouping of low-level parameters (Process Variables) as exposed through Epics Channel Access. At first, this may not seem very interesting. However, Epics Instruments allows you to

The Epics Instruments application allows you to group these components into a logical group – an Instrument. Once an Instrument has been defined, you can then save and restore settings for this Instrument. Furthermore, these settings are automatically saved in a single, portable file for later use.

Epics Channel Access gives a simple and robust interface to its lowest common unit – the Process Variable or PV. The Epics control system also provides sophisticated ways to express and manipulate complex devices, both physical and virtual. Creating such devices and defining their behavior is generally done by well-trained programmers. The application here uses a much simpler approach that can expose some categories of “Settings” that may need to be changed en masse, and returned to at a later time.

As defined here, An Epics Instrument is simply a named collection of Epics Process Variables (PVs). The PVs do not need to be physically related to one another nor be associated with a single Epics Record or Device. Rather, an Instrument is defined at the level of the Epics Channel Access client, allowing a station scientist or engineer to use their own grouping of PVs as an abstract “Instrument”. A simple example would be a pair of motors that work together to move some device.

In addition to a name and a set of PVs, an Instrument has a set of “named Positions”. At any point, the current values of an Instruments PVs can be saved as its Position. And, of course, the named Positions can then be restored simply by selecting the Position.

## 7.2 Downloading and Installation

Many Epics Applications are available as GUI Application. For Windows, you can download and install from a binary installers below.

For Linux and Mac OS X, building and installing from source are currently the principle options (a Mac OS X App will be available soon).

### 7.2.1 Downloads

The latest source kits and installers will be at


Download Option	Location
Windows Installers	<a href="#">cars_downloads</a>
Source Kit	<a href="#">github_repo</a>

## 7.3 Area Detector Display

Epics Area Detector Display is a wxPython GUI application for viewing images from an Epics Area Detector. This application requires wxPython, pyepics, numpy, and the Python Image Library.

**To run this application, simply run AD\_Display.py at the command line::** `python AD_Display.py`

and enter the base name of the AreaDetector in the PV box. A sample display would look like this:



The fields and buttons on the top control several Area Detector settings, such as starting and stopping the acquisition.

## 7.4 Strip Chart Display

StripChart is a wxPython GUI application for viewing time traces of PVs as a strip chart. It features interactive graphics, with click-and-drag zooming, updating the plotted time range, saving figures as high-quality PNGs, and saving data to ASCII files. Stripchart is inspired somewhat by the classic Epics Stripchart application written with X/Motif, but has many differences.

### 7.4.1 Dependencies, Installation

This application needs pyepics, numpy, matplotlib, and wxPython.

It also needs the wxmplot plotting library, which can be found at <http://pypi.python.org/pypi/wxmplot/>, with development versions at <http://github.com/newville/wxmplot/> and may be installed with:

```
easy_install -U wxmplot
```

Installation of the stripchart can be done with:


```
python setup.py install
```

### 7.4.2 Running Stripchart

To run this application, simply run stripchart.py at the command line:

```
python pyepics_stripchart.py
```


and enter the base name of the PVs to follow. A sample display would look like this:



### 7.4.3 Usage

Plot details can be adjusted from the configuration form, available from the Options Menu, and shown below. From this frame, you can adjust trace colors, symbols, line width and style, symbol size and styles, axes labels, and the contents and location of a plot legend. Text for titles, axes labels, and legend can include latex strings for math/Greek characters.

From the main plot, Ctrl-C works to copy to the system clipboard, and Ctrl-P will open a print dialog.



**Plot Configuration**

Title: 2011-Sep-14 10:13:02

Y Label: Py:ao2

Y2 Label: 13IDA:DMM1Ch3\_raw.VAL

X Label: Elapsed Time (seconds)

Colors: ☐ Grid ☐ Background Show: ☒ Grid ☐ Legend ☒ Legend Frame

Legend Location: upper right on plot Text Size: 12

#	Label	Color	Line Style	Thickness	Symbol	Size
1	Py:ao2	<span style="color: blue;">■</span>	solid	3	no symbol	8
2	13IDA:DMM1Ch3	<span style="color: red;">■</span>	solid	3	no symbol	8

OK

## 7.5 Using Epics Instruments

Epics Instruments is a GUI application (using wxPython) that lets any user:

- Organize PVs into Instruments: a named collection of PVs
- Manage Instruments with modern Notebook-tab interface.
- Save Positions for any Instrument by name.
- Restore Positions for any Instrument by name.
- Remember Settings for all definitions into a single file that can be loaded later.
- Multiple Users can be using multiple instrument files at any one time.

It was originally written to replace and organize the multitude of similar MEDM screens that appear at many workstations using Epics.

### 7.5.1 Dependencies, Installation

This application needs pyepics, numpy, wxPython, and SQLAlchemy (available at <http://www.sqlalchemy.org/>).

To install, use:

```
python setup.py install
```

from the Instruments folder. A windows installer may be available.


## 7.5.2 Getting Started

To run Epics Instruments, use:

```
python pyepics_instruments.py
```

or click on the icon.

A small window to select an Epics Instrument File, like this




If this is your first time using the application, choose a name, and hit return to start a new Instrument File. The next time you run Epics Instruments, it should remember which files you've recently used, and present you with a drop-down list of Instrument Files. Since all the definitions, positions, and settings are saved in a single file, restoring this file will recall the earlier session of instrument definitions and saved positions.

An Epics **Instrument** is a collection of PVs. Each Instrument will also have a collection of **Positions**, which are just the locations of all the PVs in the instrument at the time the Position was saved. Like a PV, each Instrument and each Position for an Instrument has a unique name.

## 7.5.3 Defining a New Instrument

To define a new Instrument, select **Create New Instrument** from the Instruments Menu. A screen will appear in which you can name the instrument and the PVs that belong to the Instrument.

If you add a few PVs and click OK, the PVs will connect, and you will see a screen something like this





### 7.5.4 Editing an Existing Instrument

**Edit Instrument Sample Stage**

Instrument Name:

Current PVs:	Display Type:	Remove?:
13XRM:m1.VAL	<input type="text" value="motor"/>	<input type="text" value="No"/>
13XRM:m2.VAL	<input type="text" value="motor"/>	<input type="text" value="No"/>
13XRM:m3.VAL	<input type="text" value="motor"/>	<input type="text" value="No"/>

New PVs:	Display Type	Remove?
<input type="text"/>	<input type="text"/>	<input type="text" value="No"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="No"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="No"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="No"/>
<input type="text"/>	<input type="text"/>	<input type="text" value="No"/>

### 7.5.5 The Instrument File

All the information for definitions of your Instruments and their Positions are saved in a single file – the Instruments file, with a default extension of ‘.ein’ (Epics INstruments). You can use many different Instrument Files for different domains of use.

The Instrument File is an SQLite database file, and can be browsed and manipulated with external tools. Of course, this can be a very efficient way of corrupting the data, so do this with caution. A further note of caution is to avoid having a single Instrument file open by multiple applications – this can also cause corruption. The Instrument files can be moved around and copied without problems.

## 7.6 Using Epics Motor Setup

Epics Motor Setup is a fairly simple GUI application (using wxPython) for setting up Epics Motors. A full configuration screen is shown for each motor, using a Notebook display:

Epics Motor Setup

File

13IDE:m1 13IDE:m2 13IDE:m3 13IDE:m4 13IDE:m5 ×

Label  units  Precision

**Drive**

	User	Dial	Raw	
High Limit	<input type="text" value="42.5055"/>	<input type="text" value="41.0000"/>		
Readback	<input type="text" value="38.0849"/>	<input type="text" value="36.5794"/>	<input type="text" value="414756"/>	<input type="button" value="Stop"/>
Move	<input type="text" value="38.0849"/>	<input type="text" value="36.5794"/>	<input type="text" value="414756"/>	<input type="button" value="Pause"/>
Low Limit	<input type="text" value="-73.4945"/>	<input type="text" value="-75.0000"/>	<input type="button" value="Enable"/>	<input type="button" value="Move"/>
Tweak	<input type="button" value="←"/> <input type="text" value="10.0000"/> <input type="button" value="→"/>	<input type="button" value="Disable"/>	<input type="button" value="Go"/>	

**Calibration**

Mode:   Freeze Offset:

Direction:   Offset Value:

**Dynamics**

	Normal	Backlash
Max Speed	<input type="text" value="14.1731"/>	
Speed	<input type="text" value="0.5000"/>	<input type="text" value="0.5000"/>
Base Speed	<input type="text" value="0.0050"/>	
Accel (s)	<input type="text" value="0.5000"/>	<input type="text" value="0.2000"/>
Backlash Distance		<input type="text" value="0.0000"/>
Move Fraction		<input type="text" value="1.0000"/>

**Resolution, Readback, and Retries**

Motor Res	<input type="text" value="0.0001"/>	Encoder Res	<input type="text" value="0.0001"/>
Steps / Rev	<input type="text" value="400"/>	Units / Rev	<input type="text" value="0.0353"/>
Readback Res	<input type="text" value="0.0000"/>	Readback Delay (s)	<input type="text" value="0.0000"/>
Retry Deadband	<input type="text" value="0.0001"/>	Max Retries	<input type="text" value="0"/>
Use Encoder	<input type="button" value="No"/>	Use Readback	<input type="button" value="No"/>
Use NTM	<input type="button" value="YES"/>	NTM Factor	<input type="text" value="2"/>

The real advantage of this Application is in two simple features:

1. being able to write a **motors.template** file for all the open motors.
2. being able to read and save motor settings for future use to a database. By default, a local SQLite database is used, but users at GSECARS can also use a GSE-wide database.

To save the **motors.template** information, simply type Ctrl-T to copy the full template paragraph to the system clipboard, then copy into the appropriate file. A simple example is:

```
file "$(CARS)/CARSAApp/Db/motor.db"
{
pattern
{P, M, DTYP, C, S, DESC, EGU, DIR, VELO, VBAS, ACCL, BDST,
BVEL, BACC, S REV, UREV, PREC, DHLM, DLLM}
# VAL=1.999519, OFF=0.799999, NTM=1
{13IDE:, m1, "OMS VME58", 0, -1, "Upstream Y", mm, Pos, 0.200000, 0.005000, 1.00 0000, 0.000000, 0
# VAL=2.001424, OFF=-1.000026, NTM=1
{13IDE:, m2, "OMS VME58", 0, -1, "Inboard Y", mm, Pos, 0.200000, 0.005000, 1.000 000, 0.000000, 0
# VAL=1.999519, OFF=1.069979, NTM=1
{13IDE:, m3, "OMS VME58", 0, -1, "Outboard Y", mm, Pos, 0.200000, 0.005000, 1.000000, 0.000000, 0.20
# VAL=11.035400, OFF=1.505489, NTM=1
{13IDE:, m4, "OMS VME58", 0, -1, "Upstream X", mm, Pos, 0.500000, 0.005000, 0.500000, 0.000000, 0.50
# VAL=38.084885, OFF=1.505479, NTM=1
{13IDE:, m5, "OMS VME58", 0, -1, "Downstream X", mm, Pos, 0.500000, 0.005000, 0.500000, 0.000000, 0.5
}
```

A few things to note here are:

1. You may need to change “\$(CARS)/CARSAApp/Db/motor.db” to point to the correct location of the motor.db file.
2. Though the motor “device type” and “card” are filled out, the “slot” is not available as a PV, and so is not filled out.
3. The VAL (User Value), OFF (Dial Offset) and NTM field are saved for each motor as a comment.

## 7.7 XRF Collector

XRF Collector provides a simple Epics interface control of an x-ray fluorescence detector. The main point of showing it here is as an example of using a python process to interact with a custom Epics database to provide a customized way to acquire data.

## 7.8 Ion Chamber

This application, like the XRF Collector, provides a simple Epics interface control of a non-trivial device. Again, the main point of showing it here is as an example of using a python process to interact with a custom Epics database to provide a customized way to acquire data. While the standard Epics solution would be to write a state-notation-language program that runs in an IOC, the approach here minimizes Epics coding to writing a simple database, and putting all the logic and control into a python script that runs as a long-running process along-side an IOC process.